



DENTAL INSTRUMENTS CATALOGUE



Lower Roots 33

DT-51-033

Tooth Ext Forceps



Lower Roots 01

DT-51-033-01

Tooth Ext Forceps



Lower Roots 03

DT-51-033-03

Tooth Ext Forceps



Lower Roots 06

DT-51-033-06

Tooth Ext Forceps



Incisors And Bicuspid 34N

DT-51-034

Tooth Ext Forceps

Fig -34N,



Incisors And Bicuspid 35N

DT-51-035

Tooth Ext Forceps

Fig -35N,



Tooth Ext Forceps

DT-51-036

Tooth Ext Forceps



Lower Roots 46

DT-51-046

Tooth Ext Forceps



Lower Roots 46

DT-51-046-01

Tooth Ext Forceps



Lower Roots 46

DT-51-046-02

Tooth Ext Forceps



Very Fine Roots 49

DT-51-049

Tooth Ext Forceps



Roots 51

DT-51-051

Tooth Ext Forceps



Roots 51A

DT-51-051-01

Tooth Ext Forceps



Roots 51 C

DT-51-051-03

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-057

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-057-1

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-057-02

Tooth Ext Forceps



Third Molars 67

DT-51-067

Tooth Ext Forceps



Third Molars 67 A

DT-51-067-01

Tooth Ext Forceps



Third Molars Deep Griping 67 N

DT-51-067-02

Tooth Ext Forceps



Lower Molars 73

DT-51-073

Tooth Ext Forceps

Fig -77,



Lower Roots 74

DT-51-074

Tooth Ext Forceps



Incisors And Roots 74N

DT-51-074-06

Tooth Ext Forceps

Fig -74N,



Incisors And Roots ,Slender Jaw74XN

DT-51-074-07

Tooth Ext Forceps

Fig -74XN,



Incisors And Bicuspids

DT-51-075

Tooth Ext Forceps

Fig -75,



Upper Roots 76

DT-51-076

Tooth Ext Forceps



Upper Roots 76N

DT-51-076-07

Tooth Ext Forceps



Third Molars 79

DT-51-079


Tooth Ext Forceps



Third Molars 79A

DT-51-079-01

Tooth Ext Forceps



Third Molars ,Deep Gripping 79N

DT-51-079-02

Tooth Ext Forceps



**Molars With Carious Or Broken Caps
86A**

DT-51-086-01

Tooth Ext Forceps



**Molars With Carious Or Broken Caps
86B**

DT-51-086-02

Tooth Ext Forceps



**Molars With Carious Or Broken Caps
86C**

DT-51-086-03

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-097

Tooth Ext Forceps



Very Fine Roots 133

DT-51-133

Tooth Ext Forceps



Roots Fragments, Deep Gripping 146

DT-51-146

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-151

Tooth Ext Forceps



Tooth Ext Forceps

DT-51-153

Tooth Ext Forceps



Roots, With Parallel Beak 233

DT-51-233

Tooth Ext Forceps




Incisors And Bicuspid 1

DT-51-301

Tooth Ext Forceps

Fig -1,



Incisors And Bicuspid 2

DT-51-302

Tooth Ext Forceps

Fig -2,



Bicuspid 7

DT-51-307

Tooth Ext Forceps

Fig -7,





Bicuspid 13

DT-51-313

Tooth Ext Forceps




Fig -13,



Molars 17

DT-51-317


Tooth Ext Forceps



Molars 18

DT-51-318

Tooth Ext Forceps



Molars 22

DT-51-322

Tooth Ext Forceps



Roots 33

DT-51-330


Tooth Ext Forceps



Roots 33

DT-51-333


Tooth Ext Forceps



Roots 51

DT-51-351

Tooth Ext Forceps



Roots, Very Fine Beak 51A

DT-51-351-01

Tooth Ext Forceps



Third Molars 67 A

DT-51-367-01

Tooth Ext Forceps



Molars 73

DT-51-373

Tooth Ext Forceps



Incisors And Roots 74

DT-51-374

Tooth Ext Forceps

Fig -74,




Incisors And Roots , With Fine Beak

74N

DT-51-374-01

Tooth Ext Forceps

Fig -74N,



Third Molars 79

DT-51-379

Tooth Ext Forceps



Incisors And Cuspids 1

DT-51-401

Tooth Ext Forceps

Fig -1,



Bicuspids 7

DT-51-407

Tooth Ext Forceps

Fig -7,




Incisors And Cuspids 13

DT-51-413

Tooth Ext Forceps

Fig -13,



Molars , Right 17

DT-51-417

Tooth Ext Forceps



Molars , Left 18

DT-51-418

Tooth Ext Forceps



Molars 22

DT-51-422

Tooth Ext Forceps



Roots 33

DT-51-433


Tooth Ext Forceps



Very Fine Roots 46

DT-51-446

Tooth Ext Forceps



Very Fine Roots 49

DT-51-449

Tooth Ext Forceps



Roots 51

DT-51-451

Tooth Ext Forceps



Roots , With Fine Beak 51A

DT-51-451-01

Tooth Ext Forceps



Martin Ergo 67A

DT-51-467-01

Tooth Ext Forceps



Martin Ergo 74N

DT-51-474

Tooth Ext Forceps



Martin Ergo 79

DT-51-479

Tooth Ext Forceps



Martin Ergo ,Very Fine Roots 97

DT-51-497



Bicuspid 13S

DT-51-613-09

Tooth Ext Forceps

Fig -13s,



Molars 22s

DT-51-622-09

Tooth Ext Forceps



Incisors And Roots 30S

DT-51-630-09

Tooth Ext Forceps

Fig -30S,



Roots 33S

DT-51-633-09

Tooth Ext Forceps



Incisors And Bicuspids 37

DT-51-637

Tooth Ext Forceps

Fig -37,



Incisors And Cuspids 38

DT-51-638

Tooth Ext Forceps

Fig -38,



Molars 39

DT-51-639

Tooth Ext Forceps



Molars 39A

DT-51-639-01

Tooth Ext Forceps



Molars , Left 39L

DT-51-639-05

Tooth Ext Forceps



Molars , Right 39R

DT-51-639-08

Tooth Ext Forceps



Molars 40

DT-51-640

Tooth Ext Forceps



Roots 51S

DT-51-651-09

Tooth Ext Forceps



Roots 76S

DT-51-676-09

Tooth Ext Forceps



Molars 3

DT-51-703

Tooth Ext Forceps



Incisors 5

DT-51-705

Tooth Ext Forceps

Fig -5,



Molars 6

DT-51-706

Tooth Ext Forceps



Roots 7

DT-51-707

Tooth Ext Forceps



Incisors And Cuspids 137

DT-51-737

Tooth Ext Forceps

Fig -137,



Bicuspid 139

DT-51-739

Tooth Ext Forceps

Fig -139,



Klein , Roots 51S

DT-51-751-09

Tooth Ext Forceps

Fig -51S,



Molars 17SK

DT-51-817-09

Tooth Ext Forceps



Molars 23SK

DT-51-823-09

Tooth Ext Forceps



Teeth And Roots ,Universal 150SK

DT-51-850-09

Tooth Ext Forceps



Teeth And Roots ,Universal 151Sk

DT-51-851-09

Tooth Ext Forceps




Incisors And Cuspids 1

DT-51-001

Tooth Ext Forceps

Fig -1,



Incisors And Bicuspid 2

DT-51-002

Tooth Ext Forceps

Fig -2,




Lower Incisors And Canines 4

DT-004

Tooth Ext Forceps

Fig -4,



Bicuspid 7

DT-51-007

Tooth Ext Forceps

Fig -7,




Lower Bicuspid, Universal 13

DT-51-013

Tooth Ext Forceps

Fig -13,



Molars, Right 17

DT-51-017

Tooth Ext Forceps



Molars 17SK

DT-51-817-09

Tooth Ext Forceps



Molars 23SK

DT-51-823-09

Tooth Ext Forceps



Teeth And Roots ,Universal 150SK

DT-51-850-09

Tooth Ext Forceps



Teeth And Roots ,Universal 151Sk

DT-51-851-09

Tooth Ext Forceps



Incisors And Bicuspid 1

DT-52-001

Tooth Ext Forceps Amr
18.5cm,



Henahan Incisors

DT-52-001-01

Tooth Ext Forceps Amr



Molars And Third Molars 10S

DT-52-010-09

Tooth Ext Forceps Amr
Fig -10s,



Molars 15

DT-52-015

Tooth Ext Forceps Amr



Molars 16

DT-52-016

Tooth Ext Forceps Amr
Fig -16,



Molars , Slender Jaw

DT-52-016-09

Tooth Ext Forceps Amr
Fig -16,



Molars 17

DT-52-017

Tooth Ext Forceps Amr

Fig -17,




Harris Molars , Left 18L

DT-52-018-05

Tooth Ext Forceps Amr

Fig -18L,



Harris Molars , Right 18R

DT-52-018-08

Tooth Ext Forceps Amr

Fig -18R,



Molars

DT-52-023

Tooth Ext Forceps Amr

Fig -23,



Parmlly Molars And Bicuspid 32

DT-52-032

Tooth Ext Forceps Amr

Fig -32,




Parmlly Molars And Bicuspid , Roots

32A

DT-52-032-01

Tooth Ext Forceps Amr

Fig -32L,



Molars , Left 53L

DT-52-053-05

Tooth Ext Forceps Amr
18.5cm,



Molars , Right 53R

DT-52-053-08

Tooth Ext Forceps Amr
Fig -53R,



Teeth And Roots,Universal Also For Children 62

DT-

Tooth Ext Forceps Amr
Fig -62,



Incisors And Roots 65

DT-52-065

Tooth Ext Forceps Amr
Fig -65,



Tomes Roots Fragments And Small Roots 69

DT-52-069

Tooth Ext Forceps Amr



Nevius Molars , Left 88L

DT-52-088-05

Tooth Ext Forceps Amr
Fig -88L,



Nevius Molars , Right 88R

DT-52-088-08

Tooth Ext Forceps Amr

Fig -88R,



**Cryer Bicuspid, Incisors And
Roots 150**

DT-52-150

Tooth Ext Forceps Amr



**Cryer Bicuspid, Incisors And
Roots 150A**

DT-52-150-01

Tooth Ext Forceps Amr

Fig -150A,



Bicuspid And Incisors , Roots 150X

DT-52-150-02

Tooth Ext Forceps Amr

Fig -150X,



**Bicuspid And Incisors, Roots
150AS**

DT-52-150-08

Tooth Ext Forceps Amr

fig -150 A,



**Cryer Bicuspid And Incisors ,
Roots 150S**

DT-52-150-09

Tooth Ext Forceps Amr

Fig -150S,



**Cryer Bicuspid , Incisors And Roots
151**

DT-52-151

Tooth Ext Forceps Amr

Fig -151,



**Cryer Bicuspid , Incisors And Roots
151A**

DT-52-151-01

Tooth Ext Forceps Amr

Fig -151A,



**Cryer Bicuspid , Incisors And Roots
151X**

DT-52-151-02

Tooth Ext Forceps Amr

Fig 151X,



Bicuspid And Incisors , Roots 151AS

DT-52-151-08

Tooth Ext Forceps Amr

Fig -151 AS,



**Cryer Teeth And Roots , Universal
151S**

DT-52-151-09

Tooth Ext Forceps Amr

Fig -151S,



Third Molars 222

DT-52-222

Tooth Ext Forceps Amr

Fig -222,



Root Fragments And Universal 501

DT-52-501

Tooth Ext Forceps Amr

Fig -502,



Witzel Root Fragments And Universal 502

DT-52-502

Tooth Ext Forceps Amr



Seldin 1L

DT-53-001-05


Root Elevators



Seldin 1R

DT-53-001-08


Root Elevator



Seldin 2

DT-53-002-00


Root Elevators



Seldin

DT-53-003-00

Root Elevators



Cryer

DT-53-003-01

Root Elevators



Seldin 4L

DT-53-004-05

Root Elevators



Seldin 4R

DT-53-004-08

Root Elevators



Crane Pick 6

DT-53-006-00

Root Elevators



Miller 1

DT-53-008-00

Root Elevators



Miller 2

DT-53-009-00

Root Elevators



Miller-Apexo

DT-53-010-00

Root Elevators



Miller-Apexo

DT-53-011-00

Root Elevators



Cryer
DT-53-012-00
Root Elevators




Cryer
DT-53-027-00
Root Elevators



Cryer 28
DT-028-00
Root Elevators



Seldin 34
DT-53-034-00
Root Elevators




Seldin 34S
DT-53-034-09
Root Elevators



Cryer 27S
DT-53-039-00
Root Elevators



Cryer 28S
DT-53-040-00
Root Elevators



Flohr
DT-53-046-00
Root Elevators




Barry 50
DT-53-050-00
Root Elevators
Fig -50,



Barry 51
DT-53-051-00
Root Elevators
Fig -51,



Lindo Levian 1
DT-53-200-01
Root Elevators




Lindo Levian 3
DT-53-200-03



Howard-Apexo # 1

DT-53-204-00


Root Elevators



Root Elevators


DT-53-204-01

Root Elevators



Root Elevators

DT-53-205-00



Root Elevators

DT-53-205-01

Root Elevators




Bein 1

DT-53-220-01

Root Elevator


Fig -1,




Bein 2

DT-53-220-02


Root Elevator



Bein 4
DT-53-220-04




Bein 1A
DT-53-221-01




Bein 2A
DT-53-221-02

Fig -1A,



Bein 3A
DT-53-221-03




Bein 4A
DT-53-221-04




Apical 1
DT-53-222-01

Fig -1,



Apical 2
DT-53-222-02

Fig -2,




Apical 3
DT-53-222-03

Fig -3,



Apical 4
DT-53-222-04




Apical 1A
DT-53-223-01

Fig -1A,




Apical 2A
DT-53-223-02

Fig -2,



Apical 3A
DT-53-223-03



Apical 4A
DT-53-223-04

Fig -4A,



Apical 301
DT-53-301-00

Fig -301,



Apical 302
DT-53-302-00


Root Elevators



Apical-Flohr
DT-53-303-00


Root Elevators

Fig -303,



Apical 304
DT-53-304-00

Root Elevators




Apical 46R
DT-53-311-00



Root Elevators





Coleman
DT-53-405-01
Root Elevators




Coleman 2L
DT-53-405-02
Root Elevators




Coleman 3R
DT-53-405-03
Root Elevators



Coupland Bone Chisels 1
DT-53-415-01
Root Elevators



Coupland Bone Chisels 2
DT-53-415-02
Root Elevators



Coupland Bone Chisels 3
DT-53-415-03
Root Elevators



Friedman 1
DT-53-505-01
Root Elevators




Friedman 2L
DT-53-505-02
Root Elevators



Friedman 2R
DT-53-505-03
Root Elevators




Hyllin 1
DT-53-508-01
Root Elevators



Hyllin 2
DT-53-508-02
Root Elevators




Hyllin 3L
DT-53-508-03
Root Elevators




Hylin 4R
DT-53-508-04
Root Elevators




Bein 0
DT-53-511-00
Root Elevators
Fig -0,




Bein 1
DT-53-511-01
Root Elevators
Fig -1,



Bein 2
DT-53-511-02
Root Elevators
Fig -2,



Bein 1
DT-53-512-01
Root Elevators
Fig -1,



Bein 2
DT-53-512-02
Root Elevators




Bein 0

DT-53-513-00

Root Elevators

Fig -0,




Bein 1

DT-53-513-01

Root Elevators

Fig -1,



Bein 2

DT-53-513-02

Root Elevators




Bein 0A

DT-53-514-00

Root Elevators

Fig -0A,




Bein 1A

DT-53-514-01

Root Elevators

Fig -1A,



Bein 2A

DT-53-514-02

Root Elevators




Bein 0A

DT-53-515-00

Root Elevators

Fig -0A,




Bein 1A

DT-53-515-01

Root Elevators

Fig -1A,



Bein 2A

DT-53-515-02

Root Elevators



W.F. Barry 320L

DT-53-533-20

Root Elevators



W.F. Barry 321R

DT-53-533-21

Root Elevators



Barry 322

DT-53-533-22

Root Elevators

Fig -322,



Barry 323

DT-53-533-23

Root Elevators

Fig -323,



W.F. Barry 370R

DT-53-533-70

Root Elevators



W.F. Barry 371R

DT-53-533-71

Root Elevators



Winter 11L

DT-53-551-05

Root Elevators



Winter 11R

DT-53-551-08

Root Elevators



Winter 13L

DT-53-553-05

Root Elevators



Winter 13R

DT-53-553-08

Root Elevators



Winter 14L

DT-53-554-05

Root Elevators



Winter 14R

DT-53-554-08

Root Elevators



Lecluse 1

DT-53-561-01

Root Elevators



Pott 1

DT-53-567-01

Root Elevators



Pott 2

DT-53-567-02

Root Elevators



Krallenheber 1L

DT-53-571-01

Root Elevators



Krallenheber 2R

DT-53-571-02

Root Elevators



Warwick-James 1

DT-53-576-01

Root Elevators



Warwick-James 2L

DT-53-576-02


Root Elevators



Warwick-James 3R

DT-53-576-03

Root Elevators



Chompret 1

DT-53-577-01

Root Elevators



Chompret 2

DT-53-577-02

Root Elevators



Chompret 3

DT-53-577-03

Root Elevators



Heidebrink 1

DT-53-582-01

Root Elevators



Heidebrink 2

DT-53-582-02

Root Elevators



Heidebrink 3

DT-53-582-03

Root Elevators



West 1

DT-53-584-01


Root Elevators



West 2

DT-53-584-02

Root Elevators



West 3

DT-53-584-03

Root Elevators



Bone Rongeur

DT-54-101

Bone Holding Forceps
7", 185mm,



Bone Rongeur

DT-54-102

Bone Holding Forceps
7", 185mm,



Bone Rongeur

DT-54-103

Bone Holding Forceps
8", 200mm,



Bone Rongeur

DT-54-104

Bone Holding Forceps
8", 215mm,



Bone Rongeur

DT-54-105

Bone Rongeur
145mm, 5",



Friedman-Micro 6"

DT-54-106

Bone Rongeurs
5", 150mm,



Friedman-Micro 5"

DT-54-107

Bone Rongeurs
145mm, 5",



Bone Rongeur

DT-54-108

Bone Rongeur
5", 140mm,



LUER 6"

DT-54-109

Bone Rongeurs
150mm, 6",



LUER 6"

DT-54-110

Bone Rongeurs
150mm, 6",



2A-MEAD, SQUARE TIPS 5"

DT-54-111

Bone Rongeur
145mm, 5",



BLUMENTHAL 6"

DT-54-112

Bone Rongeur
6", 155mm,



BLUMENTHAL 45° 6"

DT-54-113

Bone Rongeur
150mm, 6",



Bone Rongeur 7"

DT-54-114

Bone Rongeur
7", 175mm,



Bone Rongeur

DT-54-115

Bone Rongeur
5", 170mm,



BEYER 7"

DT-54-116

Bone Rongeur
7", 180mm,



BANE-HEARTMANN 7"

DT-54-118

Bone Rongeur
175mm, 7",



BEYER DOUBLE-ACTION 7"

DT-54-119

Bone Rongeur
7", 180mm,



ZUFAL-JANSEN 7"

DT-54-120

Bone Rongeur
7", 180mm,



Bone Rongeur

DT-54-121

Bone Rongeur
130mm, 5.5",



Bone Rongeur

DT-54-122

Bone Rongeur
130mm, 5.5",



Cleveland 5"

DT-54-123

Bone Cutting Forceps
5", 130mm,



Bone Rongeur

DT-54-124

Bone Cutting Forceps
5", 160mm,



Fig. 1
Polokoff

Polokoff

DT-55-101

Bone Files

Fig -1,



Fig. 33

Bone Files

DT-55-102

Bone Files

Fig -33,




Fig. 64
Miller

Miller

DT-55-103

Bone Files

Fig -64,



Fig. 45
Miller

Miller

DT-55-104

Bone Files

Fig -45,



Fig. 10

Bone Files

DT-55-105

Bone Files

Fig -10,



Fig. 1
Wahl

Wahl

DT-55-106

Bone Files

Fig -1,



Fig. 2
Wahl

Wahl

DT-55-107

Bone Files

Fig -2,



Fig. 12CA

Bone Files

DT-55-108

Bone Files

Fig -12CA,




Fig. 12A

Bone Files

DT-55-109

Bone Files

Fig -12A,




Fig. 12
Howard

Howard

DT-55-110

Bone Files

Fig -12,



Fig 1
Hemingway

Hemingway

DT-56-101

Bone Curettes

Fig -1,



Fig 2
Hemingway

Hemingway

DT-56-102

Hemingway

Fig -2,



Fig 3
Hemingway

Hemingway

DT-56-103

Bone Curettes

Fig -3,



Fig 4
Hemingway

Hemingway

DT-56-104

Bone Curettes

Fig -4,



Fig. 86
Lucas

Lucas

DT-56-105

Bone Curettes

Fig -86,



Fig. 87
Lucas

Lucas

DT-56-107

Bone Curettes

Fig -87,



Fig. 88
Lucas

Lucas

DT-56-108

Bone Curettes

Fig -88,



Fig. 1
Hemingway

Hemingway

DT-56-109

Bone Curettes

Fig -1,



Fig. 2
Hemingway

Hemingway

DT-56-110

Bone Curettes

Fig -2,



Fig. 3
Hemingway

Hemingway

DT-56-111

Bone Curettes



Fig. 4
Hemingway

Hemingway

DT-56-112

Bone Curettes

Fig -4,



Fig. 85
Lucas

Lucas

DT-56-113

Bone Curettes

Fig -85,



Fig. 86
Lucas

Lucas

DT-56-114

Bone Currettes

Fig -86,



Fig. 87
Lucas

Lucas

DT-56-115

Bone Currettes

Fig -87,



Fig. 88
Lucas

Lucas

DT-56-116

Bone Currettes

Fig -83,



Cotton & Dressing Tweezers

DT-57-101

Tooth Tweezers
5", 130mm,



College Tweezer

DT-57-102

Cotton & Dressing Tweezers
5", 130mm,



College Tweezer Smooth Plan Handle

DT-57-103

Cotton & Dressing Tweezers
150mm, 6",



Meriam Tweezers

DT-57-104

Cotton & Dressing Tweezers
6", 160mm,



College Tweezer

DT-57-105

Cotton & Dressing Tweezers
6", 160mm,



Meriam Tweezers

DT-57-106

Cotton & Dressing Tweezers
6", 160mm,



Meriam Tweezers

DT-57-107

Cotton & Dressing Tweezers
6", 160mm,



College Tweezer

DT-57-108

Cotton & Dressing Tweezers
150mm, 6",



College Tweezer

DT-57-109

Cotton & Dressing Tweezers
150mm, 6",



Self Locking Tweezer Serr.

DT-57-110

Cotton & Dressing Tweezers
6", 160mm,



College Tweezer

DT-57-111

Cotton & Dressing Tweezers
6", 160mm,



College Tweezer

DT-57-112

Cotton & Dressing Tweezers
150mm, 6",



Perry Tweezer

DT-57-113

Cotton & Dressing Tweezers
130mm,



Cotton & Dressing Tweezers

DT-57-114

Cotton & Dressing Tweezers
6", 160mm,



Flagg Tweezer Smooth

DT-57-115

Cotton & Dressing Tweezers
150mm, 6",



Flagg Tweezer Serrated

DT-57-116

Cotton & Dressing Tweezers
5", 160mm,



Cotton & Dressing Tweezers

DT-57-117

Cotton & Dressing Tweezers
5", 150mm,



Cotton & Dressing Tweezers

DT-57-118

Tooth Tweezers
6", 160mm,



Cotton & Dressing Tweezers

DT-57-119

Tooth Tweezers
6", 160mm,



Self Locking Tweezer Smooth

DT-57-120

Cotton & Dressing Tweezers
150mm, 6",



Suture Forceps

DT-57-121

Cotton & Dressing Tweezers
150mm, 6",



Rubber Dam Punch Ivory

DT-58-101

Rubber-Dam Instruments

6", 160mm,



Anisworth

DT-58-102

Rubber-Dam Instruments

6", 160mm,



Brewer

DT-58-103

Rubber-Dam Instruments

6", 170mm,



Palmer

DT-58-104

Rubber-Dam Instruments

7", 175mm,



Crile-Murray

DT-60-101

Needle Holders

5", 130mm,



Mayo-Hegar

DT-60-102

Needle Holders

7", 160mm,



Crile-Murray

DT-60-103

Needle Holders

150mm, 6",



Mayo-Hegar

DT-60-104

Needle Holders

150mm, 6",



Mayo-Hegar

DT-60-105

Needle Holders

7", 180mm,




Mayo-Hegar

DT-60-106

Needle Holders

8", 200mm,




Olsen-Hegar

DT-60-107

Needle Holders

5", 140mm,



Gillies

DT-60-108

Needle Holders

7", 150mm,



Boynton

DT-60-109

Needle Holders

5", 125 mm,



Mathieu

DT-60-110

Needle Holders

7", 170mm,



Mathieu

DT-60-111

Needle Holders

8", 200mm,



Crile

DT-60-112

Needle Holders

150mm, 6",



Toennis
DT-60-113
Needle Holders



Mayo-Hegar
DT-60-114
Needle Holders
7", 165mm,



Mayo-Hegar
DT-60-115
Needle Holders
8", 210mm,



Mayo-Hegar
DT-60-116
Needle Holders
165mm,



Mayo-Hegar
DT-60-117
Needle Holders
6", 165mm,



Debakey
DT-60-118
Needle Holders
5", 160mm,



Debakey
DT-60-119
Needle Holders
7",



Ryder
DT-60-120
Needle Holders



Ryder
DT-75-121
Needle Holders



Ryder 7 CARB-N-SERT Jaws
DT-60-122
Needle Holders



Bozemann 5
DT-60-123
Needle Holders
5", 135mm,



Olsen-Hegar CARB-N-SERT Jaws
DT-60-124
Needle Holders



Olsen-Hegar CARB-N-SERT Jaws

DT-60-125

Needle Holders



Olsen-Hegar CARB-N-SERT Jaws

DT-60-126

Needle Holders



Olsen-Hegar CARB-N-SERT Jaws

DT-60-127

Needle Holders



Needle Holders

DT-60-128

Needle Holders

7", 185mm,



Needle Holders

DT-60-129

Needle Holders

7", 185mm,



Ryder 8 CARB-N-SERT Jaws

DT-60-130

Needle Holders



Needle Holders

DT-60-131

Needle Holders

8", 200mm,



Mathieu 5

DT-60-132

Needle Holders

5", 140mm,



Mathieu 6

DT-60-133

Needle Holders

5", 170mm,



Mathieu 8

DT-60-134

Needle Holders

8", 200mm,



Mathieu 5

DT-60-135

Needle Holders

5", 140mm,



Mathieu 6

DT-60-136

Needle Holders

6", 170mm,



LISTER Bandage Scissors

DT-61-101

GUM Scissors

Extra fine



LISTER Bandage Scissors

DT-61-102

Gum Scissors

Extra fine



LISTER Bandage Scissors

DT-61-103

Gum Scissors

Extra fine



KNOWLES Bandage Scissors

DT-61-104

Gum Scissors

Straight



HEATH Suture Scissors

DT-61-105

Gum Scissors Serrated blade

6-1/8" (15.6cm),



SPENCER Stitch Scissors

DT-61-106

Gum Scissors

Delicate



SHORTBENT Stitch Scissors

DT-61-107

Gum Scissors Curved delicate



LITTAUER JR. Stitch Scissors

DT-61-108

Gum Scissors



NORTHBENT Stitch Scissors

DT-61-109

Gum Scissors Curved light pattern



LITTAUER Stitch Scissors

DT-61-110

Gum Scissors Standard pattern



LONG Oral Surgery Stitch Scissors

DT-61-111

Gum Scissors



MILTEX Stitch Scissors

DT-61-112

Gum Scissors Angled delicate



REEH Stitch Scissors

DT-61-113

Gum Scissors

Sharp tips small hook on 1 blade ribbon type



MILTEX Stitch Scissors

DT-61-114

Gum Scissors Curved sharp pointed tips
ribbon type



MILTEX Stitch Scissors

DT-61-115

Gum Scissors Curved sharp pointed tips
ribbon type



VANTAGE Crown Scissors

DT-61-116


Gum Scissors Collar curved



VANTAGE Collar

DT-61-117

Gum Scissors Crown Scissors straight



MILTEX Wire Cutting Scissors

DT-61-118

Gum Scissors Angled to side one serrated
blade



MILTEX Wire Cutting Scissors

DT-61-119

Gum Scissors Angled to side, one serrated blade, Carb-N-Sert inserts



Standard Pattern Operating Scissors

DT-61-120

Gum Scissors Straight sharp-sharp points



Standard Pattern Operating Scissors

DT-61-121

Gum Scissors Straight sharp-sharp points



Standard Pattern Operating Scissors

DT-61-122

Gum Scissors Straight sharp-blunt points



Standard Pattern Operating Scissors

DT-61-123

Gum Scissors Straight sharp-blunt points



Standard Pattern Operating Scissors

DT-61-124

Gum Scissors Straight sharp-blunt points



Standard Pattern Operating Scissors

DT-61-125

Gum Scissors Curved sharp/sharp points



Standard Pattern Operating Scissors

DT-61-126

Gum Scissors Curved sharp/sharp points



Standard Pattern Operating Scissors

DT-61-127

Gum Scissors Curved



Standard Pattern Operating Scissors

DT-61-128

Gum Scissors Curved, sharp-blunt points



MAYO Dissecting Scissors

DT-61-129

Gum Scissors Straight, standard beveled blades



MAYO Dissecting Scissors

DT-61-130

Gum Scissors Straight, standard beveled blades



MAYO Dissecting Scissors

DT-61-131

Gum Scissors Curved, standard beveled blades



MAYO Dissecting Scissors

DT-61-132

Gum Scissors Curved, standard beveled blades



Delicate Pattern METZENBAUM Scissors

DT-61-133

Gum Scissors Straight, blunt points
5-1/2" (14cm),



Delicate Pattern METZENBAUM Scissors

DT-61-134

Gum Scissors Straight, sharp points
5-1/2" (14cm),



Standard Pattern METZENBAUM Scissors

DT-61-135

Gum Scissors
Curved, blunt points



Standard Pattern METZENBAUM Scissors

DT-61-136

Gum Scissors
Curved, blunt points



**Delicate Pattern METZENBAUM
Scissors**

DT-61-137

Gum Scissors Curved, blunt points
5-1/2" (14cm),



**Delicate Pattern METZENBAUM
Scissors**

DT-61-138

Gum Scissors Curved, sharp points
5-1/2" (14cm),



KELLY Scissors

DT-61-139

Gum Scissors
Straight, sharp points



KELLY Scissors

DT-61-140

Gum Scissors
Curved, sharp points



KELLY Scissors

DT-61-141

Gum Scissors
Sharp points, angled on side



**KAYE Blepharoplasty And Dissecting
Scissors**

DT-61-142

Gum Scissors One serrated blade, ribbon
type
4-1/2" (10.58cm),



DEAN Dissecting Scissors

DT-61-143

Gum Scissors Blades angled on flat, one serrated blade
6-3/4" (17.5cm),



KAYE Dissecting Scissors

DT-61-144

Gum Scissors Curved, one serrated blade
4-17/32" (11.5 cm),



WAGNER Plastic Surgery Scissors

DT-61-145

Gum Scissors
Straight, sharp-sharp points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-146

Gum Scissors

Straight, sharp-sharp points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-147

Gum Scissors

Straight, sharp-blunt points



WAGNER Plastic Surgery Scissors

DT-61-148

Gum Scissors

Straight, sharp-blunt points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-149

Gum Scissors

Straight, blunt-blunt points



WAGNER Plastic Surgery Scissors

DT-61-150

Gum Scissors

Straight, blunt-blunt points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-151

Gum Scissors

Curved, sharp-sharp points



WAGNER Plastic Surgery Scissors

DT-61-152

Gum Scissors

Curved, sharp-sharp points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-153

Gum Scissors

Curved, sharp-blunt points



WAGNER Plastic Surgery Scissors

DT-61-154

Gum Scissors

Curved, sharp-blunt points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-155

Gum Scissors

Curved, blunt-blunt points



WAGNER Plastic Surgery Scissors

DT-61-156

Gum Scissors

Curved, blunt-blunt points, serrated blade



WAGNER Plastic Surgery Scissors

DT-61-157

Gum Scissors

Straight, sharp-sharp points



WAGNER Plastic Surgery Scissors

DT-61-158

Gum Scissors

Straight, sharp-sharp point, serrated blade



**Delicate Pattern METZENBAUM
Scissors**

DT-61-159

Gum Scissors Curved

5" (12.7cm),



Iris Scissors

DT-61-160

Gum Scissors Straight

4" (10.2cm),



Iris Scissors

DT-61-161

Gum Scissors Curved
4" (10.5cm),



Iris Scissors

DT-61-162

Gum Scissors Straight
4-1/2" (10.58cm),



Iris Scissors

DT-61-163

Gum Scissors Straight
4-1/2" (10.58cm),



Iris Scissors

DT-61-164

Gum Scissors Curved
4-1/2" (10.58cm),



Iris Scissors

DT-61-165

Gum Scissors Curved
4-1/2" (10.58cm),



Iris Scissors

DT-61-166

Gum Scissors Angled on side
4-1/2" (10.58cm),



Iris Scissors

DT-61-161

Gum Scissors Curved
4" (10.5cm),



Iris Scissors

DT-61-162

Gum Scissors Straight
4-1/2" (10.58cm),



Iris Scissors

DT-61-163

Gum Scissors Straight
4-1/2" (10.58cm),



Iris Scissors

DT-61-164

Gum Scissors Curved
4-1/2" (10.58cm),



Iris Scissors

DT-61-165

Gum Scissors Curved
4-1/2" (10.58cm),



Iris Scissors

DT-61-166

Gum Scissors Angled on side
4-1/2" (10.58cm),



BROPHY Scissors

DT-61-167

Gum Scissors Straight, sharp points
5-1/2" (14cm),



BROPHY Scissors

DT-61-168

Gum Scissors Curved, sharp points
5-1/2" (14cm),



GOLDMAN-FOX Scissors

DT-61-169

Gum Scissors Curved fine tips, one serrated
blade
4" (10.5cm),



GOLDMAN-FOX Scissors

DT-61-170

Gum Scissors Curved fine tips, one serrated
blade
5-1/2" (14cm),



KELLY Scissors

DT-61-171

Gum Scissors Sharp points, straight, one
serrated blade
6-1/4" (15.9 cm),



KELLY Scissors

DT-61-172

Gum Scissors
Sharp points, curved, one serrated blade



KELLY Scissors

DT-61-173

Gum Scissors

Sharp points, angular, one serrated blade



Iris Operating Scissors

DT-61-174

Gum Scissors Curved sideways

4-3/5" (14cm),



GOLDMAN-FOX Scissors

DT-61-175

Gum Scissors Curved, one serrated blade,
angled shank

4-1/2" (10.58cm),



LA GRANGE Scissors

DT-61-176

Gum Scissors

Curved, one serrated blade, angled shank



BROPHY Operating Scissors

DT-61-177

Gum Scissors Curved, sharp/sharp, with one
serrated blade

5-1/2" (14cm),



LOCKLIN Operating Scissors

DT-61-178

Gum Scissors

Angled on side 25 degrees, one serrated
blade



LOCKLIN Operating Scissors

DT-61-179

Gum Scissors

Curved shank, angled on side 25 degrees,
one serrated blade



**MILTEX SuperCut GOLDMAN-FOX
Gum Scissors**

DT-61-180

Gum Scissors

Curved



HEATH Suture Scissors

DT-61-181

Gum Scissors Serrated blade

6-1/8" (15.6cm),



Iris Scissors

DT-61-182

Gum Scissors Straight, with 20 mm blades,
delicate

3-1/2" (8.9 cm),



Iris Scissors

DT-61-183

Gum Scissors Curved, with 20 mm blades,
delicate

3-1/2" (8.9 cm),



NOYES (ROSS) Iris Scissors

DT-61-184

Gum Scissors

Curved, sharp/sharp points



CASTROVIEJO Corneal Scissors

DT-61-185

Gum Scissors Curved, sharp points
4-3/4" (.5cm),



QUIMBY Gum Scissors

DT-61-186

Gum Scissors
Strongly curved and pointed, smooth blades
15 mm



JOSEPH Scissors

DT-61-187


Gum Scissors Straight
5-1/2" (14cm),



JOSEPH Scissors

DT-61-188


Gum Scissors Curved
5-1/2" (14cm),



Hollenback

DT-62-101

Amalgam Instruments



Hollenback

DT-62-107

Amalgam Instruments



Mortonson 1

DT-62-113

Amalgam Instruments



Amalgam Instruments

DT-62-114

Amalgam Instruments

Fig -2,



Amalgam Instruments

DT-62-115

Amalgam Instruments

Fig -3,



Amalgam Instruments

DT-62-116

Amalgam Instruments

Fig -4,



Fig. 3-SHO

Amalgam Instruments

DT-62-118

Amalgam Instruments

Fig -3-SHO,



Fig. SHO-A

Amalgam Instruments

DT-62-117

Amalgam Instruments

Fig -SHO-A,



Fig. 104

Amalgam Instruments

DT-62-119

Amalgam Instruments

Fig-104,



Fig. 8

Amalgam Instruments

DT-62-120

Amalgam Instruments

Fig -8,



Fig. 1

Amalgam Instruments

DT-62-121

Amalgam Instruments

Fig -1,



Fig. 2

Amalgam Instruments

DT-62-122

Amalgam Instruments

Fig -2,



Fig 3

Amalgam Instruments

DT-62-123

Amalgam Instruments

Fig -3,



Fig 4

Amalgam Instruments

DT-62-124

Amalgam Instruments

Fig -4,



Fig SHO-A

Amalgam Instruments

DT-62-125

Amalgam Instruments

Fig -SHO-A,



Fig. 3-SHO

Amalgam Instruments

DT-62-126

Amalgam Instruments

Fig -3-SHO,



Fig. 104

Amalgam Instruments

DT-62-127

Amalgam Instruments

Fig -104,



Fig. 8

Amalgam Instruments

DT-62-128

Amalgam Instruments

Fig -8,



Fig. 1

Amalgam Instruments

DT-62-129

Amalgam Instruments

Fig -1,



Fig. 2

Amalgam Instruments

DT-62-130

Amalgam Instruments

Fig -2,



Fig. 3

Amalgam Instruments

DT-62-131

Amalgam Instruments

Fig -3,



Fig. 4

Amalgam Instruments

DT-62-132

Amalgam Instruments

Fig -4,



Fig. 1

Amalgam Instruments

DT-62-133

Amalgam Instruments

Fig -1,



Fig. 2

Amalgam Instruments

DT-62-134

Amalgam Instruments

Fig -2,



Fig. 3

Amalgam Instruments

DT-62-135

Amalgam Instruments

Fig -3,



Fig. 4

Amalgam Instruments

DT-62-136

Amalgam Instruments

Fig -4,



Ivory Fig. 8

5 mm

6 mm

7 mm

6 mm


Toffemire-Universal

Toffemire-Junior

Toffemire-Universal

DT-62-137

Amalgam Instruments



Nystrom I + II

Ivory Fig. 1

Ivory Fig. 12

Amalgam Instruments

DT-62-142

Amalgam Instruments

Levory Fig -1, Levory Fig -2,



Amalgam Instruments

DT-62-143

Amalgam Pistolen

210mm,



Amalgam Instruments

DT-62-144

Amalgam Pistolen

185mm,



Amalgam Instruments

DT-62-145

Amalgam Pistolen



Amalgam Instruments

DT-62-146

Amalgam Pistolen
145mm,



Amalgam Carrier

DT-62-147

Amalgam Carrier
Fig -1, Fig -2, Fig -4, Fig -5,



Amalgam Carrier

DT-62-148

Amalgam Carrier
Fig -1, Fig -2, Fig -3,



Amalgam Carrier

DT-62-149

Amalgam Carrier
Fig -1, Fig -2, Fig -3,



Heister

DT-63-141

Mouth Gags,
5", 130mm,



MOLT PEDIATRIC

DT-63-142

Mouth Gags,
4", 100mm,



MOLT CHILD

DT-63-143

Mouth Gags,
120 mm, 4",



MOLT ADULT

DT-63-144

Mouth Gags,
9", 140mm,



FERGUSON-ACKLAND

DT-63-145

Mouth Gags,
145mm, 5",



DENHART

DT-63-146

Mouth Gags,
5", 125 mm,



ROSER-KONIG

DT-63-148

Mouth Gags,
7", 185mm,



LANGENBECK

DT-63-149

Finger Protectors
7", 180mm,



Tongue Depressor And Cheek Retractors

DT-63-150

Tongue Depressor and Cheek Retractors



Retractor

DT-63-151

Mouth Gags,



HARTMANN

DT-63-152

Retractors
145mm, 5",



Retractor

DT-63-153

Mouth Gags,
5", 140mm,



SCHMID

DT-63-154

Mouth Gags,



Retractor

DT-63-155

Mouth Props
5", 140mm,



Retractor

DT-63-156

Retractors
6", 160mm,



VESTIBULUM RETRACTOR

DT-63-157

Retractors



STERNBERG

DT-63-158

Retractors



ROUX

DT-63-159

Retractors
23cm 36 x 30,



ROUX
DT-63-160
Retractors



ROUX
DT-63-161
Retractors



MIDDELDORPF
DT-63-162
Retractors
200mm,



MIDDELDORPF
DT-63-163
Retractors
8", 210mm,



MIDDELDORPF
DT-63-164
Cheek Retractors
9", 230mm,



Gauze Case
DT-63-165
Gauze Case
9", 26x30mm, 230mm,



Retractor

DT-63-166

Retractors
5", 165mm,



Retractor

DT-63-167

Retractor
6", 165mm,



Retractor

DT-63-168

Retractor
6", 165mm,



Retractor

DT-63-169

Retractor
6", 165mm,



Retractor

DT-63-170

Retractor
6", 165mm,



Retractor

DT-63-171

Retractor
6", 165mm,



Retractor

DT-63-172


Retractors
6", 165mm,



Retractor

DT-63-173

Retractors
6", 165mm,



Retractor

DT-63-174

Retractors
6", 160mm,



Retractor

DT-63-175

Retractors
6", 160mm,



LANGENBECK-GREEN

DT-63-176

Retractors



Fig. 72
1962/1965
Kirkland

Kirkland
DT-64-101

Periodontia Instruments



Fig. 112
Orban

Orban
DT-64-102

Periodontia Instruments



Fig. 112
Orban

Orban
DT-64-104

Periodontia Instruments



Fig. 7
Goldman-Fox

Goldman-Fox
DT-64-105

Periodontia Instruments



Fig. 8
Goldman-Fox

Goldman-Fox
DT-64-106

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox
DT-64-107

Periodontia Instruments



Fig. 1-2
Kirkland

Kirkland
DT-64-103

Periodontia Instruments



Fig. 11
Goldman-Fox

Goldman-Fox
DT-64-108

Periodontia Instruments



Fig. 3/4
Buck

Buck
DT-64-109

Periodontia Instruments

Fig -3/4,



Fig. 5/6
Buck

Buck
DT-64-110

Periodontia Instruments

Fig -5/6,



Fig. 7
Goldman-Fox

Goldman-Fox
DT-64-111

Periodontia Instruments



Fig. 8
Goldman-Fox

Goldman-Fox
DT-64-112

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox

DT-64-113

Periodontia Instruments



Fig. 11
Goldman-Fox

Goldman-Fox

DT-64-114

Periodontia Instruments



Fig. 3/4
Buck

Buck

DT-64-115

Periodontia Instruments

Fig -3/4,



Fig. 5/6
Buck

Buck

DT-64-116

Periodontia Instruments

Fig -5/6,



Fig. 1/2
Sanders

Sanders

DT-64-117

Periodontia Instruments



Fig. 4/5
Sanders

Sanders

DT-64-118

Periodontia Instruments



Fig. 18/20

OCHSENBEIN CHISEL NO 4

DT-64-119

Periodontia Instruments



Fig. 5/6
Back

Crane-Kaplan

DT-64-120

Periodontia Instruments

Fig -5/6,



Fig. 1/2
Sander

Sanders

DT-64-121

Periodontia Instruments



Fig. 4/5
Sander

Sanders

DT-64-122

Periodontia Instruments



Fig. 18/20

OCHSENBEIN CHISEL NO 4

DT-64-123

Periodontal Instruments



Fig. CK3
Crane-Kaplan

Crane-Kaplan

DT-64-124

Periodontia Instruments

Fig -CK3,



Fig. HU 1/2
Curette

Curette

DT-64-125

Periodontia Instruments

Fig -HU 1/2,



Fig. 2R/2L
Columbia

Columbia

DT-64-126

Periodontia Instruments

Fig -2R/2L,



Fig. 4R/4L
Columbia

Columbia

DT-64-127

Periodontia Instruments

Fig -4A/4L,



Fig. 13/14
Columbia

Columbia

DT-64-128

Periodontia Instruments

Fig -13/14,



Fig. 411/412

Periodontia Instruments

DT-64-129

Periodontia Instruments



Fig. 413/414

Periodontia Instruments

DT-64-130

Periodontia Instruments



Zerling-Chiles

Zerling-Chiles

DT-64-131

Periodontia Instruments



Fig. HU 1/2
Curette

Curette

DT-64-132

Periodontia Instruments

Fig -HU 1/2,



Fig. 2R/2L
Columbia

Columbia

DT-64-133

Periodontia Instruments

Fig -2R/2L,



Fig. 4R/4L
Columbia

Columbia

DT-64-134

Periodontia Instruments

Fig -4A/4L,



Fig. 13/14
Columbia

Columbia

DT-64-135

Periodontia Instruments

Fig -13/14,



Fig. 1/2
Langer

Langer

DT-64-136

Periodontia Instruments



Fig. 3/4
Langer

Langer

DT-64-137

Periodontia Instruments



Fig. 1/2
Barnhart

Barnhart

DT-64-138

Periodontia Instruments

Fig -1/2,



Fig. 5/6
Barnhart

Barnhart

DT-64-139

Periodontia Instruments

Fig -5/6,



Fig. ME-1
(10x11)

Periodontia Instruments

DT-64-140

Periodontia Instruments



Fig. 3/4
ME-1 (23/24)

Periodontia Instruments

DT-64-141

Periodontia Instruments



Fig. 1/2
Langer

Langer

DT-64-142

Periodontia Instruments



Fig. 3/4
Langer

Langer

DT-64-143

Periodontia Instruments



Fig. 1/2
Barnhart

Barnhart

DT-64-144

Periodontia Instruments

Fig -1/2,



Fig. 5/6
Barnhart

Barnhart

DT-64-145

Periodontia Instruments

Fig -5/6,



Fig. ME-1
[10/11]

Periodontia Instruments

DT-64-146

Periodontia Instruments



Fig. 3/4
Rue [23/24]

Periodontia Instruments

DT-64-147

Periodontia Instruments



Fig. 1/2
Gracey

Gracey

DT-64-148

Periodontia Instruments



Fig. 3/4
Gracey

Gracey

DT-64-149

Periodontia Instruments



Fig. 5/6
Gracey

Gracey

DT-64-150

Periodontia Instruments



Fig. 7/8
Gracey

Gracey

DT-64-151

Periodontia Instruments



Fig. 9/10
Gracey

Gracey

DT-64-152

Periodontia Instruments



Fig. 11/12
Gracey

Gracey

DT-64-153

Periodontia Instruments



Fig. 13/14
Gracey

Gracey

DT-64-154

Periodontia Instruments



Fig. 1/2
Gracey

Gracey

DT-64-155

Periodontia Instruments



Fig. 3/4
Gracey

Gracey

DT-64-156

Periodontia Instruments



Fig. 5/6
Gracey

Gracey

DT-64-157

Periodontia Instruments



Fig. 7/8
Gracey

Gracey

DT-64-158

Periodontia Instruments



Fig. 9/10
Gracey

Gracey

DT-64-159

Periodontia Instruments



Fig. 11/12
Gracey

Gracey

DT-64-160

Periodontia Instruments



Fig. 13/14
Gracey

Gracey

DT-64-161

Periodontia Instruments



Fig. 11/14
Gracey

Gracey

DT-64-162

Periodontia Instruments



Fig. 12/13
Gracey

Gracey

DT-64-163

Periodontia Instruments



Fig. 1
De Marco

Ce Marco

DT-64-164

Periodontia Instruments

Fig -1,



Fig. 2
De Marco

De Marco

DT-64-165

Periodontia Instruments

Fig -1,



Fig. 13/14
Anterior

Anterior

DT-64-166

Periodontia Instruments

Fig -13/14,



Fig. 17/18
Posterior

Posterior
DT-64-167

Periodontia Instruments



Fig. 11/14
Gracey

Gracey
DT-64-168

Periodontia Instruments



Fig. 12/13
Gracey

Gracey
DT-64-169

Periodontia Instruments



Fig. 1
De Marco

De Marco
DT-64-170

Periodontia Instruments

Fig -1,



Fig 2
De Marco

De Marco
DT-64-171

Periodontia Instruments

Fig -2,



Fig. 13/14
Anterior

Anterior
DT-64-172

Periodontia Instruments

Fig -13/14,



Fig. 17/18
Poserior

Poserior

DT-64-173

Periodontia Instruments



Fig. 1
Goldman-Fox

Goldman-Fox

DT-64-174

Periodontia Instruments



Fig. 2
Goldman-Fox

Goldman-Fox

DT-64-175

Periodontia Instruments



Fig. 3
Goldman-Fox

Goldman-Fox

DT-64-176

Periodontia Instruments



Fig. 4
Goldman-Fox

Goldman-Fox

DT-64-177

Periodontia Instruments



Fig. 5
Goldman-Fox

Goldman-Fox

DT-64-178

Periodontia Instruments



Fig. 6
Goldman-Fox

Goldman-Fox

DT-64-179

Periodontia Instruments



Fig. 21
Goldman-Fox

Goldman-Fox

DT-64-180

Periodontia Instruments



Fig. 1
Goldman-Fox

Goldman-Fox

DT-64-181

Periodontia Instruments



Fig. 2
Goldman-Fox

Goldman-Fox

DT-64-182

Periodontia Instruments



Fig. 3
Goldman-Fox

Goldman-Fox

DT-64-183

Goldman-Fox



Fig. 4
Goldman-Fox

Goldman-Fox

DT-64-184

Periodontia Instruments



FIG. 5
Goldman-Fox

Goldman-Fox

DT-64-185

Periodontia Instruments



Fig. 6
Goldman-Fox

Goldman-Fox

DT-64-186

Periodontia Instruments



Fig. 21
Goldman-Fox

Goldman-Fox

DT-64-187

Periodontia Instruments



Fig. 11/12
Mc Call

Mc Call

DT-64-188

Periodontia Instruments



Fig. 11A/12A
Mc Call

Mc Call

DT-64-189

Periodontia Instruments



Fig. 13/14
Mc Call

Mc Call

DT-64-190

Periodontia Instruments



Fig. 13S/14S
Mc Call

Mc Call

DT-64-191

Periodontia Instruments



Fig. 17S/18S
Mc Call

Mc Call

DT-64-192

Periodontia Instruments



Fig. 19/20
Mc Call

Mc Call

DT-64-193

Periodontia Instruments



Fig. 11/12
Mc Call

Mc Call

DT-64-194

Periodontia Instruments



Fig. 11A/12A
Mc Call

Mc Call

DT-64-195

Periodontia Instruments



Fig. 13/14
Mc Call

Mc Call

DT-64-196

Periodontia Instruments



Fig. 13S/14S
Mc Call

Mc Call

DT-64-197

Periodontia Instruments



Fig. 1/S/18S
Mc Call

Mc Call

DT-64-198

Periodontia Instruments



Fig. 1B/20
Mc Call

Mc Call

DT-64-199

Periodontia Instruments



Fig. 2/3
Taylor

Taylor

DT-64-200

Periodontia Instruments



Fig. 204S

Periodontia Instruments

DT-64-202

Periodontia Instruments



Fig. 204S

Periodontia Instruments

DT-64-203

Periodontia Instruments



Fig. 6
Crane-Kaplan

Crane-Kaplan

DT-64-204

Periodontia Instruments

Fig -6,



Fig. 77/0
Younger-Good

Younger-Good

DT-64-205

Periodontia Instruments



Fig. 273
Taylor

Taylor

DT-64-206

Periodontia Instruments



Fig. 204

Periodontia Instruments

DT-64-207

Periodontia Instruments



Fig. 204S

Periodontia Instruments

DT-64-208

Periodontia Instruments



Fig. 204S

Periodontia Instruments

DT-64-209

Periodontia Instruments



Fig. 6
Crane-Kaplan

Crane-Kaplan

DT-64-210

Periodontia Instruments

Fig -6,



Fig. 7/8
Younger-Good

Younger-Good

DT-64-211

Periodontia Instruments



Fig. 3/7
Hirschfeld

Hirschfeld

DT-64-212

Periodontia Instruments



Fig. 28/33
Jaquette

Jaquette

DT-64-213

Periodontia Instruments



Fig. 31/32
Jaquette

Jaquette

DT-64-214

Periodontia Instruments



Fig. 31/32
Jaquette

Jaquette

DT-64-215

Periodontia Instruments



Fig. 3/7
Hirschfeld

Hirschfeld

DT-64-216

Periodontia Instruments



Fig. 30/33
Jaquette

Jaquette

DT-64-217

Periodontia Instruments



Fig. 31/32
Jaquette

Jaquette

DT-64-218

Periodontia Instruments



Fig. 34/35
Jaquette

Jaquette

DT-64-219

Periodontia Instruments



Fig. 37/38
Jaquette

Jaquette

DT-64-220

Periodontia Instruments



Fig. 1
Ochsenbein

Ochsenbein

DT-64-221

Periodontia Instruments



Fig. 2
Ochsenbein

Ochsenbein

DT-64-222

Periodontia Instruments



Fig. 1
Gardner

Gardner

DT-64-223

Periodontia Instruments

Fig -1,



Fig. 2
Gardner

Gardner

DT-64-224

Periodontia Instruments



Fig. 5
Gardner

Gardner

DT-64-225

Periodontia Instruments



Fig. 5R
Gardner

Gardner

DT-64-226

Periodontia Instruments



Fig. 6
Gardner

Gardner

DT-64-227

Periodontia Instruments



Fig. 1W
Wakefield

Wakef Old

DT-64-228

Periodontia Instruments



Fig. 2W
Wakefield

Wakef Old

DT-64-229

Periodontia Instruments



Fig. 3
Buckley

Buckley

DT-64-230

Periodontia Instruments

Fig -3,



Fig. M3
Molt

Molt

DT-64-231

Periodontia Instruments



Fig. PR3
Prichard

Prichard

DT-64-232

Periodontia Instruments



Fig.
Hopkins

Hopkins

DT-64-233

Periodontia Instruments



Fig 2
Hopkins

Hopkins

DT-64-234

Periodontia Instruments



Fig -5
Frerr

Frerr

DT-64-235

Periodontia Instruments

Fig -5,



Fig -16
Frerr

Frerr

DT-64-236

Periodontia Instruments

Fig -16,



Fig.3
Mead

Mead

DT-64-237

Periodontia Instruments



Fig. 23
Seldin Retractor

Seldin Retractor

DT-64-238

Periodontia Instruments



Williger
c/n. 17

Williger

DT-64-239

Periodontia Instruments



Williger
cm. 16

Williger
DT-64-240

Periodontia Instruments



Williger
cm. 16

Williger
DT-64-241

Periodontia Instruments



Williger
cm. 13

Williger
DT-64-242

Periodontia Instruments



Fig. 1S/2S
Sugarman

Sugarman
DT-64-243

Periodontia Instruments



Fig. 3S/4S
Sugarman

Sugarman
DT-64-244

Periodontia Instruments



Fig. 1/2
TG Chsel

TG Chsel
DT-64-245

Periodontia Instruments



Fig. 1/2
TG Chisel

TG Chsel
DT-64-246

Periodontia Instruments



Fig. 1
TG Chisel

TG Chsel
DT-64-247

Periodontia Instruments



Fig. 2
TG Chisel

TG Chsel
DT-64-248

Periodontia Instruments



Fig. 1
(15K)
Kirkland

Kirkland
DT-64-249

Periodontia Instruments



Fig. 2
(16K)
Kirkland

Kirkland
DT-64-250

Periodontia Instruments



Fig. 1
Orban

Orban
DT-64-251

Periodontia Instruments




Fig. 2
Orban

Orban

DT-64-252

Periodontia Instruments




Westcott

DT-65-101

Filling instruments

Fig -1,




Westcott

DT-65-110

Filling instruments

Fig -1,




Filling Instruments

DT-65-119

Filling instruments

Fig -1, Fig -2, Fig -4, Fig -5, Fig -3,



Filling Instruments

DT-65-125

Filling instruments

Fig -1, Fig -2, Fig -3,



No. 22 Cement Spatula

DT-66-101

No. 22 Cement Spatula



No. 24A Cement Spatula, Rigid

DT-66-102

No. 24A Cement Spatula, Rigid



No. 24 Cement Spatula, Flexible

DT-66-103

No. 24 Cement Spatula, Flexible



No. 313 Cement Spatula,

DT-66-104

No. 313 Cement Spatula,



No. 5 Cement Spatula

DT-66-105


No. 5 Cement Spatula



No. 7 Cement Spatula

DT-66-106

No. 7 Cement Spatula




Single Ended / Double Ended Explorers

DT-67-101

Single Ended / Double Ended Explorers

Fig -UT 1/2,



Single Ended / Double Ended Explorers

DT-67-112

Single Ended / Double Ended Explorers

Fig -UT 1/2,



SIMPLE STEM PLANE 4

DT-68-113

MOUTH MIRROR



SIMPLE STEM PLANE 5

DT-68-114

MOUTH MIRROR



SIMPLE STEM PLANE 4

DT-68-115

MOUTH MIRROR



SIMPLE STEM PLANE 4

DT-68-116

MOUTH MIRROR



SIMPLE STEM PLANE 4

DT-68-117

MOUTH MIRROR



SECO / CONE-SOCKET

DT-68-118

Mirror Handles

5", 125 mm,



OCTOGONAL CONE SOCKET

DT-68-119

Mirror Handles
120 mm, 4",



MARTIN ERGO

DT-68-120

Mirror Handles
5", 135mm,



Mirror Handles

DT-68-121

Mirror Handles
5", 135mm,



Mirror Handles

DT-68-122

Mirror Handle
5", 135mm,



Mirror Handles

DT-68-123

Mirror Handles
5", 135mm,



Mirror Handles

DT-68-124

Mirror Handles
5", 135mm,



Napkin Holder

DT-68-131

Napkin Holder

ø 50 mm,



Laryngeal Mirrors

DT-68-129

MIRROR HANDLES

7 mm 30°, ø 4,



Cotton Wool Container

DT-68-130

Cotton Wool Container

ø 50 mm,



Laryngeal Mirrors

DT-68-128

MIRROR HANDLES

7 mm 30°, ø 4,



MIRROR HANDLE CONE SOCKET

DT-53-721-01

MIRROR HANDLES

14cm/5,



MIRROR HANDLE CONE SOCKET

DT-53-721-02

MIRROR HANDLE CONE SOCKET

14cm/5,



MIRROR HANDLE CONE SOCKET

DT-53-723-01

MIRROR HANDLES

13cm/5,



MIRROR HANDLE

DT-53-728-02

MIRROR HANDLES

12cm,



MARTIN-GRIP

DT-53-746-02

12cm,



MIRROR HANDLE

DT-53-755-02

MIRROR HANDLES

14cm,



MIRROR HANDLE

DT-53-757-08

MIRROR HANDLES

14cm/5,



MIRROR

DT-53-702-22



MIRROR
DT-53-703-22




RHODIUM MIRROR
DT-53-704-10
MIRROR



RHODIUM MIRROR
DT-53-704-18
MIRROR



CONE SOCKET
DT-53-706-22
MIRROR



CONE SOCKET
DT-53-707-22
MIRROR



MICRO-MOUTH MIRROR
DT-53-710-03


Fig -1,



MICRO-MOUTH MIRROR

DT-53-710-05

Fig -1,



RHODIUM MICRO MIRROR

DT-53-711-03

Fig -1,



HILGER

DT-53-715-00

18cm /7,



HILGER

DT-53-715-01

CHEEK RETRACTORS



Fig. 1
Spreader

Root Spreader Fig 1

DT-69-101

Endodontic Instruments

Fig -1,



Fig. 2
Spreader

Root Spreader Fig 2

DT-69-102

Endodontic Instruments

Fig -2,



Fig. 2
Spreader

Root Spreader Fig 3

DT-69-103

Endodontic Instruments

Fig -3,



Fig. 4
Spreader

Root Spreader Fig 4

DT-69-104

Endodontic Instruments

Fig -4,



Fig. 5
Spreader

Root Spreader Fig 5

DT-69-105

Endodontic Instruments

Fig -5,



Fig. 6
Spreader

Root Spreader Fig 6

DT-69-106

Endodontic Instruments

Fig -6,



Fig. 7
Spreader

Root Spreader Fig 7
DT-69-107

Endodontic Instruments
Fig -7,



Fig. 1
Luks

Luks Fig 1
DT-69-108

Endodontic Instruments
Fig -1,



Fig. 2

LUKS Fig 2
DT-69-109

Endodontic Instruments
Fig -2,



Fig. 3

LUKS Fig 3
DT-69-110

Endodontic Instruments
Fig -3,



Fig. 4

LUKS Fig 4
DT-69-111

Endodontic Instruments
Fig -4,



Fig. 1

Vertical Condenser Fig 1
DT-69-112

Endodontic Instruments
Fig -1,



Fig. 2

Vertical Condenser Fig 2

DT-69-113

Endodontic Instruments

Fig -2,



Fig. 3

Vertical Condenser Fig 3

DT-69-114

Endodontic Instruments

Fig -3,



Fig. 6/1

Root Canal Explorer Fig 6/1

DT-69-116

Endodontic Instruments

Fig -6/1,



Fig. 6/17

Root Canal Explorer Fig 6/17

DT-69-117

Endodontic Instruments

Fig -6/17,



Fig. 6/23

Root Canal Explorer Fig 6/23

DT-69-118

Endodontic Instruments

Fig -6/23,



Fig. 31L

Root Canal Excavators Fig 31L

DT-69-119

Endodontic Instruments

Fig -31,



Fig. 32L

Root Canal Excavators Fig 32L

DT-69-120

Endodontic Instruments

Fig -32L,



Fig. 33L

Root Canal Excavators Fig 33L

DT-69-121

Endodontic Instruments

Fig -33L,



Fig. 6/1

Root Canal Explorers Fig 6/1

DT-69-122

Endodontic Instruments

Fig -6/1,



Fig. 6/17

Root Canal Explorers Fig 6/17

DT-69-123

Endodontic Instruments

Fig -6/17,



Fig. 6/23

Root Canal Explorers Fig 6/23

DT-69-124

Endodontic Instruments

Fig -6/23,



Fig. 31L

Root Canal Excavators Fig 31L

DT-69-125

Endodontic Instruments

Fig -31L,



Fig. 32L

Root Canal Excavators Fig 32L

DT-69-126

Endodontic Instruments

Fig -32L,



Fig. 33L

Root Canal Excavators Fig 33L

DT-69-127

Endodontic Instruments

Fig -33L,



Heidemann

DT-70-101

Separating Spatulas

Fig -1,



Heidemann

DT-70-102

Separating Spatulas

Fig -2,



Heidemann

DT-70-104

Separating Spatulas

Fig -4,



Heidemann

DT-70-105

Separating Spatulas

Fig -5,



Heidemann

DT-70-106

Separating Spatulas

Fig -6,



Heidemann

DT-70-107

Separating Spatulas

Fig -7,



Heidemann
DT-70-108

Separating Spatulas
Fig -8,



Heidemann
DT-70-109

Separating Spatulas
Fig -9,



Heidemann
DT-70-110

Separating Spatulas
Fig -1,



Heidemann
DT-70-111

Separating Spatulas
Fig -2,



Heidemann
DT-70-112

Separating Spatulas
Fig -3,



Fig. 0

Excavators Fig 0

DT-71-101

Excavators



Fig. 1

Excavators Fig 1

DT-71-102

Excavators

Fig -1,



Fig. 2

Excavators Fig 2

DT-71-103

Excavators

Fig -2,



Fig. 3

Excavators Fig 3

DT-71-104

Excavators

Fig -3,



Fig. 4

Excavators Fig 4

DT-71-105

Excavators

Fig -4,



Fig. 5

Excavators Fig 5

DT-71-106

Excavators

Fig -5,



Fig 6

Excavators Fig 6

DT-71-107

Excavators

Fig -6,



Fig. 0

Excavators Fig 0

DT-71-108

Excavators



Fig. 1

Excavators Fig 1

DT-71-109

Excavators

Fig -1,



Fig. 2

Excavators Fig 2

DT-71-110

Excavators

Fig -2,



Fig. 3

Excavators Fig 3

DT-71-111

Excavators

Fig -3,



Fig. 4

Excavators Fig 4

DT-71-112

Excavators

Fig -4,



Fig. 5

Excavators Fig 5

DT-71-113
Excavators
Fig -5,



Fig. 6

Excavators Fig 6

DT-71-114
Excavators
Fig -6,



Fig. 6

Excavators Fig 6

DT-71-115
Excavators
Fig -6,



Fig. 14

Excavators Fig 14

DT-71-116
Excavators
Fig -14,



Fig. 17

Excavators Fig 17

DT-71-117
Excavators
Fig -17,



Fig. 17W

Excavators Fig 17W

DT-71-118
Excavators
Fig -17W,



Fig. 18

Excavators Fig 18

DT-71-119
Excavators
Fig -18,



Fig. 18W
Oval spoon

Oval Spoon Fig 18W

DT-71-120
Excavators



Fig. 19

Excavators Fig 19

DT-71-121
Excavators
Fig -19,



Fig. 19W
Oval spoon

Oval Spoon Fig 19W

DT-71-122
Excavators



Fig. 6

Excavators Fig 6

DT-71-123
Excavators
Fig -6,



Fig. 14

Excavators Fig 14

DT-71-124
Excavators
Fig -14,



Fig. 17

Excavators Fig 17

DT-71-125
Excavators
Fig -17,



Fig. 17W
Oval spoon

Oval Spoon Fig 17W

DT-71-126
71-126



Fig. 18W
Oval spoon

Oval Spoon 18W

DT-71-127
Excavators



Fig. 19

Excavators Fig 19

DT-71-129
Excavators
Fig -19,



Fig. 19W
Oval spoon

Oval Spoon Fig 19W

DT-71-130
Excavators



Fig. 38/39

Excavators Fig 38/39

DT-71-131
Excavators
Fig -38/39,



Fig. 63/64

Excavators Fig 63/64

DT-71-133
Excavators
Fig 63/64,



Fig. 65/66

Excavators Fig 65/66

DT-71-134
Excavators



Fig. 31L

Excavators Fig 31L

DT-71-135
Excavators
Fig -31L,



Fig. 32L

Excavators Fig 32L

DT-71-136
Excavators
Fig -32L,



Cleoid / Discoid

DT-71-137
No. 3/6 Excavator, double end



Cleoid / Discoid

DT-71-138
No. 89/92 Carver - Discoid " Cleoid, octagonal handle



Fig. U15
Towner

Towner Fig U15

DT-72-104

Scalers

Fig -U15,



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33

DT-72-105

Scalers

Fig- U15/33,



Fig. U15/30
Towner Jaquette

Towner Jaquette Fig- 15/30

DT-72-106

Scalers

Fig- 15/30,



Fig H5/33

Fig- H5/33

DT-72-107

Scalers

Fig -H5/33,



Fig. H6/H7

Fig- H6/H7

DT-72-108

Scalers

Fig- H6/H7,



Fig. 1/2

Fig- 1/2

DT-72-109

Scalers

Fig -1/2,



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33

DT-72-110

Scalers

Fig- U15/33,



Fig. U15/30
Towner Jaquette

Towner Jaquette Fig- U15/30

DT-72-111

Scalers

Fig- U15/30,




Fig. 1

Fig- 1

DT-72-112

Scalers

Fig-1,



Fig. 33

Scalers

DT-72-117

Scalers

Fig -33, Fig -31, Fig -32, Fig -34, Fig -35,



Fig. H5
Hygienist

Hygienist Fig- H5

DT-72-122

Scalers

Fig -H5,



P.K. THOMAS CARVERS FIG. 1

DT-73-101

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.1



P.K. THOMAS CARVERS FIG. 2

DT-73-102

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.2



P.K. THOMAS CARVERS FIG. 3

DT-73-103

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.3



P.K. THOMAS CARVERS FIG. 4

DT-73-104

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.4



P.K. THOMAS CARVERS FIG. 5

DT-73-105

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.5



P.K.THOMAS INSTRUMENTS

DT-73-106

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.5 SILVER



P.K.THOMAS INSTRUMENTS

DT-73-107

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.4 RED



P.K.THOMAS INSTRUMENTS

DT-73-108

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.3 BLUE



P.K.THOMAS INSTRUMENTS

DT-73-109

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.2 GOLD



P.K.THOMAS INSTRUMENTS

DT-73-110

P.K.THOMAS WAX AND MODELLING CARVERS

FIG.1 GREEN



FAHNENSTOCK Wax Carver \\\\" Knife, 5\\\\" , Small

DT-73-111

FAHNENSTOCK Wax Carver \\\\" Knife, 5\\\\" ,
small



FAHNENSTOCK Wax Carver \\\\" Knife, 6-3/4\\\\" , Large

DT-73-112

FAHNENSTOCK Wax Carver \\\\" Knife, 6-
3/4\\\\" , large



Impression Trays

DT-74-101

Impression Trays

66x55mm,

Impression Trays

DT-74-102

Impression Trays

62x50mm,

Impression Trays

DT-74-103

Impression Trays

64x50mm,



Impression Trays

DT-74-104

Impression Trays

65x46mm,

Impression Trays

DT-74-105

Impression Trays

69x54mm,

Impression Trays

DT-74-106

Impression Trays

65x46mm,



Impression Trays

DT-74-107

Impression Trays

70x48mm,



Impression Trays

DT-74-108

Impression Trays

70x55mm,



Impression Trays

DT-74-109

Impression Trays

59x42mm,



Impression Trays

DT-74-110

Impression Trays

55x37mm,



Impression Trays

DT-74-111

Impression Trays

61x42mm,



Impression Trays

DT-74-112

Impression Trays

44mm,



Impression Trays

DT-74-113

Impression Trays

44mm,



Impression Trays

DT-74-114

Impression Trays

44mm,



Impression Trays

DT-74-115

Impression Trays

66x55mm,



Impression Trays

DT-74-116

Impression Trays

65x46mm,



Impression Trays

DT-74-117

Impression Trays

69x54mm,



Impression Trays

DT-74-118

Impression Trays

69x54mm,



Impression Trays

DT-74-119

Impression Trays

70x48mm,



Impression Trays

DT-74-120

Impression Trays

70x55mm,



Impression Trays

DT-74-121

Impression Trays

55x37mm,



Impression Trays

DT-74-122

Impression Trays

61x42mm,



Impression Trays

DT-74-123

Impression Trays

61x42mm,



Impression Trays

DT-74-124

Impression Trays

61x42mm,



Impression Trays

DT-74-125

Impression Trays



Dental Syringes

DT-75-101

Dental Syringes
55x37mm,



Dental Syringes

DT-75-102

Dental Syringes
55x37mm,



Self Aspirating 1.8cc

DT-75-103

Dental Syringes
55x37mm,



Cylinder Cartridge Syringe

DT-75-104

Dental Syringes
55x37mm,



Dental Syringes

DT-75-106

Dental Syringes
55x37mm,



Dental Syringes

DT-75-107

Dental Syringes
55x37mm,



C-W 1.8cc

DT-75-108

Dental Syringes
55x37mm,



Dental Syringes

DT-75-109

Dental Syringes
55x37mm,



Dental Syringes

DT-75-110

Dental Syringes
55x37mm,



Dental Syringes

DT-75-111

Dental Syringes
55x37mm,



Dental Syringes

DT-75-112

Dental Syringes
55x37mm,



Dental Syringes

DT-75-113

Dental Syringes
55x37mm,



Dental Syringes

DT-75-118

Dental Syringes

Fig -1, Fig -2, Fig -4, Fig -3,



Dental Syringes

DT-75-119

Aspira-Plus

1.8 ml,



Dental Syringes

DT-75-120

Dental Syringes

1.8 ml,



Dental Syringes

DT-75-122

Dental Syringes

1.8 ml,



Dental Syringes

DT-75-124

Dental Syringes

1.8 ml,



Dental Syringes

DT-75-126

Dental Syringes

1, 8 ml,



Dental Syringes

DT-75-127

Dental Syringes
1, 8 ml,



Dental Syringe

DT-75-128

Aspira-Plus
1.8 ml,



Dental Syringe

DT-75-129

Dental Syringe



Dental Syringe

DT-

Dental Syringe



Dental Syringe

DT-75-131


Dental Syringe



Dental Syringe

DT-75-132

Dental Syringe
1 ml,




Wedeistaet

DT-76-101

Cavity Preparation

Fig -1, Fig -2, Fig -4, Fig -5, Fig -6, Fig -3,



Wedeistaet

DT-76-112

Cavity Preparation

Fig -1, Fig -2, Fig -4, Fig -5, Fig -6, Fig -3,



11 Briault D/E Explorer
DT-77-101

FIG. # 11 BRIAULT D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E Explorer
DT-77-102

FIG. # 11/12 OLD DOMINION UNIVERSITY D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E After Five Explorer
DT-77-103

FIG. # 11/12 ODU D/E AFTER FIVE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2A Pigtail D/E Explorer
DT-77-104

FIG. # 2A PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2H D/E Explorer
DT-77-105

FIG. # 2H DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2R/2L D/E Explorer
DT-77-106

FIG. # 2R/2L D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3 Pigtail D/E Explorer
DT-77-107

FIG. # 3 PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3CH Cowhorn D/E Explorer
DT-77-108

FIG. # 3CH COWHORN D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3ES D/E Explorer
DT-77-109

FIG. # 3ES DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



5 D/E Explorer
DT-77-110

FIG. # 5 D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



6/3 D/E Explorer
DT-77-111

FIG. # 6/3 DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



17 S/E Explorer
DT-77-112

FIG. # 17 SE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



23 Shepherd's Hook Explorer

DT-77-113

FIG. # 23 SHEPHERD'S HOOK S/E EXPLORER.
MADE OF SOLID STAINLESS STEEL HANDLE



3 S/E Explorer

DT-77-114

FIG. # 3 SE EXPLORER. MADE OF SOLID
STAINLESS STEEL HANDLE



3A S/E Explorer

DT-77-115

FIG. # 3A S/E EXPLORER. MADE OF SOLID
STAINLESS STEEL HANDLE



54 S/E Explorer

DT-77-116

FIG. # 54 S/E EXPLORER. MADE OF SOLID
STAINLESS STEEL HANDLE



A Tips For MORREL Crown Remover

DT-78-101

A tips for MORREL Crown Remover



Tips For MORREL Crown Remover

DT-78-102

tips for MORREL Crown Remover



C Tips For MORREL Crown Remover

DT-78-103

C tips for MORREL Crown Remover



MORREL Crown Remover, With 3 Points

DT-78-104

MORREL Crown Remover, with 3 points



WYNMAN Carbide Crown Gripper, 5-1/2", Without Ratchet, With TC Jaws

DT-78-105

WYNMAN Carbide Crown Gripper, 5-1/2", without ratchet, with TC jaws



X-Ray Film Hangers 6 Clips

DT-79-101

X-RAY FILM HANGERS 6 CLIPS MADE OF STAINLESS STEEL



X-Ray Film Hangers 10 Clips

DT-79-102

X-RAY FILM HANGERS 10 CLIPS MADE OF STAINLESS STEEL



Matrix Retainer Toffilemire
DT-80-101

MATRIX BAND RETAINER TOFFILEMIRE.
AVAILABLE IN STRAIGHT & CONTRA.
STAINLESS STEEL



Toffilemire Matrix Retainer Junior /
Child
DT-80-102

MATRIX BAND RETAINER TOFFILEMIRE
JUNIOR / CHILD. AVAILABLE IN STRAIGHT &
CONTRA. STAINLESS STEEL



Toffilemire Matrix Retainer Senior /
Adult
DT-80-103

MATRIX BAND RETAINER TOFFILEMIRE
SENIOR / CHILD. AVAILABLE IN STRAIGHT &
CONTRA. STAINLESS STEEL



Toffilemire Matrix Retainer
Universal
DT-80-104



Matrix Band Retainer Ivory 8
DT-80-105

MATRIX BAND RETAINER IVORY 8, AVAIALBE



Matrix Band Retainer Ivory Plain
DT-80-106

MATRIX RETAINER IVORY PLAIN / WITHOUT



Matrix Retainer Ivory With Shield
DT-80-107

MATRIX RETAINER IVORY WITH LATERAL SHIELD (JAW). STAINLESS STEEL



Matrix Retainer Cheek Protector
DT-80-108

MATRIX BAND RETAINER CHEEK PROTECTOR. STAINLESS STEEL



Matrix Band Retainer Nystrom Left
DT-80-109

MATRIX BAND RETAINER NYSTROM LEFT. STAINLESS STEEL



Matrix Band Retainer Nystrom Right
DT-80-110

MATRIX BAND RETAINER NYSTROM RIGHT. STAINLESS STEEL



Matrix Band Retainer Siqveland (Ash Narrow)
DT-80-111

MATRIX BAND RETAINER SIQVELAND (ASH NARROW). STAINLESS STEEL



Matrix Band Retainer Siqveland (Ash Wide)
DT-80-112

MATRIX BAND RETAINER SIQVELAND (ASH WIDE). STAINLESS STEEL



**Matrix Band Retainer Ivory Didner
(Adjustable)**

DT-80-113

Matrix Band Retainer Ivory Didner
(Adjustable)



Matrix Retainer Ivory Separator

DT-80-114

MATRIX BAND RETAINER IVORY SEPARATOR.
STAINLESS STEEL



Matrix Retainer Elliott Separator

DT-80-115

MATRIX RETAINER ELLIOTE SEPARATOR.
STAINLESS STEEL



Matrix Retainer Ivory 8 N

DT-80-116


MATRIX BAND RETAINER IVORY 8N. MADE OF
STAINLESS STEEL



Matrix Retainer Didner 8 P

DT-80-117

MATRIX BAND RETAINER IVORY 8P. MADE OF
STAINLESS STEEL



**PLASTERLESS ARTICULATOR GALETTI
(WHITE)**

DT-81-101

HIGH TECH PLASTERLESS ARTICULATOR GALETTI WITH ADJUSTABLE SCREWS AT REARS & FRONT. SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. ALL ANATOMICAL MOVEMENTS POSSIBLE



Vibratory Polished Articulator Galetti

DT-81-102

GALETTI ARTICULATOR IS SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. EASY TO USE, ALLOWING QUICK FIXATION OF MODEL. IMMEDIATE ADJUSTMENT OF THE OCCLUSAL & CORRECTS IMULATION OF CENTRAL RELATION. ALL



ANATOMICAL ARTICULATOR (ASA)

DT-81-103

ANATOMICAL ARTICULATOR WITH BEARING ROD & INCISIVE PIN, 4 MOUNTING PLATES & SUPPORTING BACK-PIN FOR 45° VIEW, OPENING 125°, INTERCONDYLAR DISTANCE 110 MM - SOLIDLY CONSTRUCTED OF METAL ALLOY CAS



UNIVERSAL ARTICULATOR

DT-81-104

ANATOMICAL ARTICULATOR, MADE OF ALUMINUM



Gnatus Pro JP30 Articulator

DT-81-105

The Gnatus Pro JP30 Articulator offers the best



ANATOMICAL ARTICULATOR

DT-81-106

ANATOMICAL ARTICULATOR WITH INCISIVE PIN, 4



KEYSTONE BLACK ARTICULATOR WIRED

DT-81-107


KEYSTONE DELUXE UNIVERSAL TYPE-A ARTICULATOR. INTERCHANGEABLE, INCISAL PIN INCLUDED, LATERAL EXCURSIONS.



KEYSTONE DELUXE MAGNET ARTICULATOR

DT-81-108

KEYSTONE DELUXE MAGNETIC DENTURE ARTICULATOR WITH PLATES, INTERCHANGEABLE, INCISAL PIN INCLUDED. MAGNETIC SPLIT CASTING ENSURES EASY & EXACT SETTING, LATERAL EXCURSIONS



PROFESSIONAL OCCLUSOR ARTICULATOR

DT-81-109

PROFESSIONAL OCCLUSOR, VIBRATORY POLISHED ARTICULATOR - ALUMINUM, BRASS, STAINLESS STEEL




Magicart-1 Articulator

DT-81-110



FULL DENTURE ARTICULATOR (ROUND)

DT-81-111



MODEL 'M' HI TOP STEPHAN ARTICULATOR

DT-81-112



Chromium Plated Plain Round

DT-81-113

Nickle Chromium Plated Plain Round Articulator



UNIVERSAL ARTICULATOR (YELLOW)

DT-81-114

UNIVERSAL ARTICULATOR OPENING 125° WITH 4 MOUNTING PLATES - CONDYLAR INCLINATION FIXED AT 30°, INTERCONDYLAR DISTANCE AT 110MM, BENNET ANGLE FIXED AT 7°, DISTANCE BETWEEN MOUNTING PLATES 85 MM.



FOSTER DENTURE ARTICULATOR

DT-81-115

THIS ARTICULATOR IS EQUIPPED WITH AN ADJUSTABLE ANTERIOR PIN STOP & POSTERIOR SCREW STOP, AND CONSTRUCTED OF A SLIGHTLY HEAVIER GAUGE BRASS GIVING THE ARTICULATOR A MORE SOLID FEEL.



FORK CHROME PLATED ARTICULATOR

DT-81-116

NICKLE CHROMIUM PLATED ARTICULATOR FORK



HINGE ARTICULATOR FULL ARCH

DT-81-117

AN ECONOMICAL & SIMPLE CHROMIUM PLATED



**QUADRANT / CROWNED BRIDGE
ARTICULATOR**

DT-81-118



UNIVERSAL ARTICULATOR (BLUE)

DT-81-119

UNIVERSAL ARTICULATOR OPENING 125° WITH 4 MOUNTING PLATES - CONDYLAR INCLINATION FIXED AT 30°, INTERCONDYLAR DISTANCE AT 110MM, BENNET ANGLE FIXED AT 7°, DISTANCE BETWEEN MOUNTING PLATES 85 MM.



FOSTER / APEX # 2 ARTICULATOR

DT-81-120

FOSTER / APEX # 2 ARTICULATOR, THE WIDTH AT THE JOINT IS 1-3/4" WITH A VERTICAL PARALLEL HEIGHT OF 2-1/8". THE SMALL SIZE OF THIS ARTICULATOR ALLOWS THE VERSATILITY TO USE FOR FULL OR PARTIAL



MORTIMER-I OCULADER ARTICULATOR

DT-81-121

PLAINLINE BALL BEARING MORTIMER-I OCULADER ARTICULATOR. ALL FILLED WITH HIGH CARBON STEEL BALLS FOR SMOOTH SETTING & ARTICULATION OF DENTURE



MORTIMER-II OCULADER ARTICULATOR

DT-81-122

PLAINLINE BALL BEARING MORTIMER-II OCULADER



OCCLUSOR, WIRE ARTICULATOR

DT-81-123

OCCLUSOR, WIRE ARTICULATOR - STAINLESS STEEL



ASA Articulator 5002

DT-81-124

ASA ARTICULATOR WITH LATERAL MOMENT. MADE



Model "A" Crescent Stephan Articulator

DT-81-125

CRESCENT STEPHAN MODEL "A" ARTICULATOR IS USED FOR PARTIALS, CENTER SPRING WITH ADJUSTABLE TENSION, SMOOTH BRASS CASTINGS WITH NICKEL PLATING, CENTER BAR IN BOTH UPPER & LOWER BOW FOR MAXIMUM



Bite Fork / Occlusal Alignment Plane

DT-81-126

BITE FORK IS A RELIABLE INSTRUMENT FOR CHECKING THE CORRECT POSITION OF THE MASTICATORY PLANE, ACHIEVING PARALLELISM OF THE OCCLUSAL PLANE WITH THE CAMPER'S LINE & THE BIPUPILLARY LINE IS FUNDAM



Bandage And Utility Scissors

DT-82-101

Bandage and Utility Scissors, 7-1/2" (19.1 cm), needle destroyer, serrated blade, autoclavable



Plaster And Compound Knife

DT-82-102

Plaster and Compound Knife



Soldering Tweezers, Curved

DT-82-103

Soldering Tweezers, curved



SWISS Jeweler Style Forceps

DT-82-104

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 3, narrow, fine, 4-3/4"



SWISS Jeweler Style Forceps

DT-82-105

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 5, Super fine, 4-3/8" (11.2



SWISS Jeweler Style Forceps

DT-82-106

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 7, fine, curved, 4-1/2"



Calipers

DT-82-107

Calipers



Calipers

DT-82-108

Calipers



Steel Ruler

DT-82-109

Steel Ruler



Forceps For Articulation Paper

DT-82-110

Forceps for articulation paper



CONTACT US

-  Street No.6 Haji Noorpura College Road Daska
-  +92 335 0044230
-  info@dentita.com
-  www.dentita.com